

SARTORIUS

Simplifying Progress



Incucyte® SX5 Training
Mathias Lucas
Field Application Scientist, France



Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion




Incucyte® Live Cell Monitoring and Analysis




Acquire




Process



Data



Incucyte® Key Applications



Cell Health & Proliferation

- Cell Counting
- Viability
- Apoptosis
- Cytotoxicity
- Tumor Spheroid
- Cell Cycle
- ATP Metabolism
- Mitochondrial Membrane Potential
- Organoid Culture QC

Cell Function

- Immune Cell Killing
- Immune Cell Activation
- Antibody Internalization
- Phagocytosis
- NETosis
- Neuronal Activity
- Live-Cell Immunocytochemistry

Cell Movement & Morphology

- Scratch Wound Migration
- Scratch Wound Invasion
- Chemotaxis
- Neurite Outgrowth
- Spheroid Invasion
- Dilution Cloning

Continuously monitor & analyze cells inside your incubator using simple mix-and-read protocols!

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion



Incucyte® SX5 Overview

3 Optical Module Options

- All modules have HD Phase (brightfield with Spheroid and Organoid scan types)

- Green/Orange/NIR




Channel	Excitation	Emission
Green	453 - 485 nm	494 - 533 nm
Orange	546 - 568 nm	576 - 639 nm
NIR	648 - 674 nm	685 - 756 nm

- Green/Red



Channel	Excitation	Emission
Green	441 - 481 nm	503 - 544 nm
Red	567 - 607 nm	622 - 704 nm



Incucyte® SX5 Overview

3 Optical Modules

- Metabolism 
- Dual Excitation, Single Emission for ATP Analysis


Channel	Excitation	Emission
Green	473 - 498 nm	NA
Orange	524 - 550 nm	565 - 591 nm
NIR	648 - 674 nm	685 - 756 nm



Incucyte® SX5 Overview


Objectives

- 4x, 10x, 20x on an automated turret
- Compatible with most tissue culture vessels :
 - Multiwell plates, dishes, flasks and slides
 - >600 vessels from >30 manufacturers
- **Software**
 - Remote access
 - Guided interface
 - Unlimited user licenses



Swapping the Optical Modules

- Detect the installed optical module in the Device menu
- Physically swap the module and then detect again to update to the new Optical Module
- Swapping modules will cancel upcoming scans



Remove the Front Cover of the Incucyte® SX5 Gantry

- Make sure that the Incucyte® SX5 is powered off
- Locate the two Cover Tabs on the bottom of the Front Cover and push outward to disengage Cover from Incucyte® SX5
- Pull Cover down and away from the Incucyte® SX5 to remove
- Place Cover to the side and in a safe location



Remove the Incucyte® SX5 Optical Module Assembly

- Identify the Optical Module and Tab on the front of the Optical Module
- Flip Tab up and turn 90 degrees counterclockwise until the detent feature disengages
- Pull the Optical Module away and out of Incucyte® SX5
- Place Optical Module in a secure location or back into the proper packaging. The Pelican™ case which contains the optional Optical Modules is fitted to hold any Optical Modules designed for the Incucyte® SX5



Remove the Incucyte® SX5 Phase Lamp

- Identify the Phase Lamp by looking for the top-right cut-out in the Front Support, where the Label on the Phase Lamp will be in-view
- Manually pull out the Drawer to the furthest extent and remove any Trays or Vessels
- Reach up and into the top of the Incucyte® SX5 towards the Label and grasp the Phase Lamp, front and back
- Pull the Phase Lamp to the right side of the Incucyte® SX5 – in the direction of the arrow
- Carefully, remove the Phase Lamp from the Incucyte® SX5
- Place Phase Lamp in a secure location or back into the proper packaging




Install the Incucyte® SX5 Phase Lamp

- Remove the new Phase Lamp from box or packaging
- With the Drawer still in the pulled-out position, guide the Phase Lamp into position
- Gently guide and push the Phase Lamp towards the left side of the Incucyte® SX5 until an audible and tactile click is noted
- The Phase Lamp will not seat into the proper position if not properly aligned



Install the Incucyte® SX5 Optical Module Assembly

- Remove new Optical Module from box or packaging
- Flip Tab up and turn 90 degrees clockwise or counterclockwise so that the Tab is vertically aligned
- Guide the Optical Module into the same position within the Incucyte® SX5 that the prior Optical Module was removed
- While gently holding the Optical Module in position, flip up and turn the tab 90 degrees clockwise until detent engages
- Flip Tab down and into the proper recess in the Optical Module
- Replace any Trays or Vessels that were removed earlier and manually push the Drawer back into the Incucyte® SX5



Replace the Front Cover of the Incucyte® SX5 Gantry

- Guide the Front Cover back into position on the Incucyte® SX5 by leading with the top and guiding upward, using the Guide Hooks for proper placement and orientation
- Once the Cover is in place and flush with the rest of the Incucyte® SX5, push the two Cover Tabs inward to engage Front Cover
- Power Controller on and allow the Incucyte® SX5 to run the warmup sequence before beginning any scanning




Incucyte® SX5 Function and Best Practices

- Incucyte® SX5 is NOT an incubator and will not regulate the environment
- Fans draw air into the gantry so conditions inside the Incucyte® are the same as the incubator but there is no active environmental control
- Operating Environmental Conditions
 - 0°C to 42°C
 - 5% to 95% Relative Humidity, Non-Condensing
- Best Practices
 - Set Incubator Temp 0.5°C below desired temp (e.g. if 37°C is desired, then set incubator to 36.5°C)
 - Check water pan and humidity settings 2-3 times per week
 - Always wear gloves and practice sterile technique

Incucyte® SX5 Maintenance


- The exterior of the Incucyte® gantry can be cleaned with 70% alcohol during routine incubator cleaning
- Vaporized multi-purpose disinfectants (e.g., Virkon) or acidic solutions can be corrosive to machinery and should NOT be used to clean the incubator
- Follow the manufacturer's recommendations for routine incubator maintenance but do NOT heat sterilize, if this is needed, remove the Incucyte® gantry from the incubator and move the controller if it is on top of the incubator
- The Incucyte® can be sterilized with vaporized hydrogen peroxide treatment
- Calibration is recommended every 6 months



Automated Image Acquisition and Auto-Focus

- Image based auto-focus
- Performs a series of wide (coarse and fine) sweeps to assess several image planes
- Finds plane of highest contrast and collects phase, green, and red image at the same plane
- Auto-focuses on every well, every plate, and every time point
- Seeding negative control wells to the upper left corner of a microwell plate can be beneficial





Pipetting Rules to Remember

- Hold the pipette vertically, when aspirating
- Immerse the tip only 2-3 mm
- Pre-rinse 3-5 times before pipetting (forward)
- Pause consistently after aspiration
- Dispense at a 30-45° angle
- Pipette against the inside wall of the receiving vessel
- Operate with smooth and consistent thumb action
- Reverse pipetting significantly reduces bubble formation



Accurately Pipette Cells and Samples: Technique is Key for Accurate Results

Tacta Mechanical Pipettes



Picus Electronic Pipettes



- Pipette with accuracy and comfort to ensure:
- Cells are evenly distributed
- Cell densities are accurate
- Volumes of treatments are accurate for dose response analysis
- Set-up high throughput assays rapidly and efficiently

Cultures should be Free of Contaminants: Technique is Key for Accurate Results

Microsart® AMP Mycoplasma Kit



MD8 Air Samplers



MYCAP™



- Check cell culture routinely for mycoplasma contamination
- Ensure cell culture environment has the lowest concentration of airborne contaminants
- Aseptically remove samples and close containers without contaminating contents

Weigh, Filter, & Concentrate: Technique is Key for Accurate Results

Sartoclear Dynamics® Lab Kits



Vivaspin® Turbo



Vivaflow®



MA37 Infrared Moisture Analyzer



Cubis® MSA



Quintix® & Secura®



- Filter cell cultures and media to ensure samples are free of unwanted debris
- Concentrate proteins, antibodies, and conditioned media from mammalian cell culture and test effects on cellular function
- Accurately weigh samples to ensure downstream assay set-up and analysis are optimal

Ultrapure Water for Multiple Applications: Technique is Key for Accurate Results



- Adapt an ultrapure water system for different applications
- Ensure water for cell culture applications is free of impurities

Live Cell Experimental Optimization

Live Cell Experimental Optimization

Condition	Things to consider
Cell Seeding	<ul style="list-style-type: none"> • Density may need to be optimized (always seed cells at least 10% confluence) • Allow cells to settle at ambient temperature for 20 minutes before placing into the incubator for optimal cell distribution • Non-adherent cells may require a coating for optimal cell distribution
Media	<ul style="list-style-type: none"> • Ensure sufficient medium per well (especially 96 and 384-well plates) • Schedule media changes if longer than 3 days • No need to use phenol red-free medium • High riboflavin (e.g. Fischer's) can give green autofluorescence
Reagents	<ul style="list-style-type: none"> • Follow protocol recommendations for optimization • Make up stock solution containing reagents to dilute treatments
Treatments	<ul style="list-style-type: none"> • Always run a positive and negative control • Run at least 3 replicates for most applications
Plate Type	<ul style="list-style-type: none"> • Choose plates optimised for the experiment, black walled are not necessary • Good optical clarity assists focusing and makes analysis simpler
Imaging	<ul style="list-style-type: none"> • Always verify data with images • Minimize scanning frequency and fluorescent channel exposure for sensitive cell types


Best Practices for Image Acquisition

Image Artefacts to Avoid	How to Get the Best Image
Empty Wells (No Cells)	Edit the scan pattern to exclude wells that do not contain cells. Ensuring >10% cell confluency in all wells is recommended (unless scanning in dilution cloning mode)
Media Droplets on the Lid	Carefully wipe away
Fingerprints/Scratches/Debris	Avoid touching tops/bottoms of vessels (even with gloved hands) and using vessels with scratches/debris
Writing on the Vessel	Label vessels on the sides or in areas where imaging will not occur
Condensation	Allow the plate to warm-up for ~20min in the incubator/ Incucyte before scanning or use a pre-warmed lid
Vessels Not Being Positioned Correctly	Take trays out of the Incucyte and place on a flat surface. Lock microplates into place with well A1 in top left corner
Bubbles	Remove bubbles with ethanol vapor (De-Bubbler / Bubble Buster). Remove the inner straw from a wash bottle, add 70% ethanol and gently squeeze vapor over the wells to clear bubbles

Scratches



Condensation



The De-Bubbler /
Bubble Buster



Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion




Live Cell Analysis Software Training

4 Steps to Image Analysis


Acquire Images

Scan Properties	
Vessel Type	24-well Corning Falcon
Scan Type	Standard
Image Channels	<input checked="" type="checkbox"/> Phase <input checked="" type="checkbox"/> Green
	Acquisition Time (ms) 300
Objective	20x
Scan Duration	2 min (estimated)
Number of Daily Scans	24


View Images



Process Images




Export Data




Incucyte® : A Guided Interface


Schedule Scans




Vessel View




Analysis Definition



Graphing and Exporting




Logging into the Incucyte® SX5



- Incucyte®: Instrument identifier assigned by your FSE at installation (e.g. IC70014 or an IP address, 192.168.128.130 for a direct connection)
- User Name and Passwords can be set up by the Admin
- Main Window opens up to guide you through the interface

Navigating the IncuCyte® Software



Data Management and Archiving

- The Incucyte® SX5 has 27.3TB of storage capacity across 4 hard drives
- Data is stored in a RAID format so is duplicated across different hard drives
- However, routine data backup is encouraged (Incucyte Manual page 151)
- Incucyte® data can be archived via the network or a directly connected hard drive
- Once archived, data can be retained on the Incucyte® SX5 or deleted
- Deleted data cannot be restored to the Incucyte® SX5 , however, reanalysis of the archived data set can be performed on the users own PC through the Incucyte® GUI
- Routine database management and removal of older experiments is recommended periodically

User Accounts

Function	Admin	Standard	Limited	Guest
Add/Remove/Edit Scheduled Scans	Any	Any	No	No
View Vessels	Any	Any	Any	Any
Create Analysis Definitions	Any	Any	Any	No
View Analysis Definitions	Any	Any	Any	Any
Edit/Delete Analysis Definitions	Any	User Owned	User Owned	No
Edit Vessel Documentation	Any	Any	Any	No
Delete Vessels	Any	User Owned	No	No
Archive Vessels	Any	Any	Any	No

Individual user accounts are recommended

Admin users can also perform software updates, load licence updates and run calibrations

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion



Getting Started – Inserting a Vessel

Device Status	LED Colour	LED State
Idle	Green	Pulsing
Active	Red	Solid
Pre-Scan (15 Seconds Before a Scan)	Yellow	Flashing
Drawer Open	Yellow	Flashing
Front Panel Removed	Red	Pulsing
Device Error	Red	Flashing



Opening the gantry unit and inserting a vessel

- Check (Software, touch display and LED-state) before opening
- Opening whilst scanning will lose that timepoint, but not the whole experiment

Acquire


Scheduling Scans

4 Steps to Image Analysis


Acquire Images

Scan Properties	
Vessel Type	24-well Corning Falcon
Scan Type	Standard
Image Channels	<input checked="" type="checkbox"/> Phase <input checked="" type="checkbox"/> Green
	Acquisition Time (ms) 300
Objective	20x
Scan Duration	2 min (estimated)
Number of Daily Scans	24


View Images



Process Images




Export Data




Incucyte®: A Guided Interface


Schedule Scans




Vessel View



Analysis Definition




Graphing and Exporting



Acquire

Scheduling Scans: Launch the Guided Interface



Launch the Guided Interface

- Click the blue + icon
- Will open a new window

Acquire

Scheduling Scans: Repeating or One Off?

Scan Repeatedly or Once?

Specify whether the vessel will be scanned repeatedly or only once.

Scan on Schedule
Select this option to add a vessel to the repeating scan schedule.

Scan Once Now
Select this option to immediately scan a single vessel just once.

Back Next

Scan Repeatedly or Once?

- Most scans will be on a repeating schedule
- Scan Once Now is useful for a quick look at a vessel, e.g., for assessing confluence

Acquire

Scheduling Scans: New or Copied Vessel?

Create or Restore Vessel

Either create a new vessel from scratch, copy an existing vessel, or restore a previously-scanned vessel.

Create Vessel

- + New**
Select this option to create a brand new vessel to scan.
- Copy Current**
Select this option to create a new vessel by copying a vessel from the current schedule.
- Copy Previous**
Select this option to create a new vessel by copying a previously scanned vessel.

Restore Vessel

- Add Scans**
Select this option to restore a previously-scanned vessel for additional scans.


Back Next

Create or Restore Vessel

- Create a new vessel from scratch
- Copy a vessel currently in the Incucyte or one scanned previously, e.g. a replicate plate
- Restore a previously scanned vessel and add scans to that same experiment

Acquire

Scheduling Scans: Types of Scan



Scan Type

- Standard scanning is used for a large number of assays including cell health assays and to select Cell-by-Cell scanning
- Application specific scan types will automatically pre-select appropriate vessels / channels / objectives

Acquire

Scheduling Scans: Channels and Objectives

Scan Settings

- Cell-by-Cell scanning, if available, can be enabled at this point
- Each vessel can be scanned using a different selection of channels
- Default acquisition times will work for most applications
- Choice of objective is also vessel specific (some imaging modes require specific objectives)

Acquire

Scheduling Scans: Choose a Vessel

IC Add Vessel

Vessel Selection

Select the type of vessel to scan.

Enter text to search...

Manufacturer	Category	Wells	Area	Catalog Numbers	Vessel Name	Tray Name
Aurora	Plate	384	N/A	ABE2-10100A, ABE2-11100A, ABE2-11101A	384-well Aurora IQ-EB Ultra Low Base	Microplates
Cellvis	Plate	12	N/A	P12-1.5H-N	12-well Cellvis	Microplates
Cellvis	Plate	96	N/A	P96-1.5H-N	96-well Cellvis	Microplates
Corning	Plate	6	N/A	3335, 3471, 3506, 3516	6-well Corning	Microplates
Corning	Plate	6	N/A	353046, 353224, 353846, 353934, 354400, 354402, 354404, 354413, 354417, 354428, 354431, 354432, 354510, 354515, 354595, 354603, 354652, 354658, 356400, 356413, 356515, 356652	6-well Corning Falcon	Microplates
Corning	Plate	12	N/A	3336, 3512, 3513	12-well Corning	Microplates
Corning	Plate	12	N/A	353043, 353225, 354470, 354500, 354501, 354502, 354503, 356470	12-well Corning Falcon	Microplates
Corning	Plate	24	N/A	3337, 3473, 3524, 3526, 3527	24-well Corning	Microplates
Corning	Plate	24	N/A	353047, 353226, 353847, 353935, 354408, 354411, 354412, 354414, 354433, 354605, 354619, 354635, 354659, 356408, 356414	24-well Corning Falcon	Microplates
Corning	Plate	48	N/A	3338, 3548	48-well Corning	Microplates
Corning	Plate	48	N/A	353078, 353230, 354505, 354506, 354507, 354508, 354509, 356505, 356509	48-well Corning Falcon	Microplates
Corning	Plate	96	N/A	3300, 3474, 3585, 3595, 3596, 3598, 3599, 3628, 3997	96-well Corning	Microplates
Corning	Plate	96	N/A	3340, 3603, 3610, 3841, 3842, 3843, 3903, 3904	96-well Corning (Blk/Wht)	Microplates
Corning	Plate	96	N/A	7007	96-well Corning Round Bottom ULA	Microplates
Corning	Plate	96	N/A	353219, 354640, 354649, 354650, 354651, 356640, 356649, 356650, 356651, 356692, 356693, 356700	96-well Corning Falcon Blk/Clr Imagi...	Microplates

Back Next

Vessel Selection

- Choose from a wide range of different tissue culture vessels
- Use the search tool to find your vessel of choice
- Some applications require specific vessels


Acquire

Scheduling Scans: Choose the Location

IC Add Vessel

Vessel Location

Specify the location in the drawer for the vessel.



Front of Instrument

Select the cutout where the vessel will be scanned.

The drawing to the left depicts the IncuCyte drawer as viewed from above.

If the schedule currently contains another vessel or tray in the selected location, a warning will appear below.

☐ Remove the following vessel(s) from the schedule:

Vessel ID	Vessel Name	Owner	First Scan
1000...	dasda	AdminUser	N/A

Back

Next

Vessel Location

- Physically load a vessel into the instrument, then schedule
- A warning will appear if trying to add a vessel to an occupied location

Acquire

Scheduling Scans: Where to Scan

IC Add Vessel

Scan Pattern

Specify the scan pattern to use for image acquisition.

Select Wells

	1	2	3	4	5	6	7	8	9	10	11	12
A	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
B	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
C	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
D	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
E	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
F	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
G	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
H	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Images per Well

Scan Duration 3 min (estimated)

[Use Sample Pattern](#)

[Back](#) [Next](#)

Scan Pattern

- Select wells containing cells, do not scan empty wells
- In a 96 well plate with the 10x objective, capture 2 images per well for adherent / homogeneous cultures
- Capture 4 images per well for non-adherent, non-homogeneous cells or with Cell-by-Cell imaging

Acquire

Scheduling Scans: Platemap & Notebook

IC Add Vessel

Vessel Notebook

Provide information about the vessel.

Name *

Cell Type

Passage

Plate Map

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

No Plate Map Specified

Notes

Back Next

Vessel Notebook

- Label every experiment clearly (make the name searchable)
- Add a Plate Map to record experimental design and automatically group replicates

Acquire

Scheduling Scans: Analysis with Acquisition

IC Add Vessel

Analysis Setup

Define an analysis to launch automatically after the vessel has been scanned. This is optional, but recommended.

Analysis Type <Defer analysis until later> ▾

Analysis Definition <Defer analysis until later>

Spectral Unmixing

Basic Analyzer

Angiogenesis

NeuroTrack

Adherent Cell-by-Cell

% G contributes to O	0 ▾
% O contributes to N	0 ▾
% G contributes to N	0 ▾
% N contributes to G	0 ▾

Analysis Notes

Back

Next

Analysis Setup

- Add analysis in parallel with image acquisition
- Only when an Analysis Definition as already been established
- Fluorescent images will need to be spectrally unmixed with previously determined values

Acquire

Scheduling Scans: Set the Scan Frequency

IC Add Vessel

Scan Schedule

Define the scan schedule for this vessel.

00 02 04 06 08 10 12 14 16 18 20 22 00

Add Scans to Schedule

- ☒ Create new schedule with scans at intervals of 3 Hours
- ☐ Add to existing schedule
- ☐ Create new schedule with advanced scheduling options
- ☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

- ☒ Scan indefinitely
- ☐ Stop scanning 1 days, 0 hours after the first scan
- ☐ Stop scanning after 11/11/2020 11:00

For advanced scheduling options, double click the timeline on the main schedule page.

Back Next

Schedule Scans

- Create new schedule, use advanced scheduling options or defer scanning
- Choose indefinite scanning or select an end time
- Fine adjustments can be made by click and drag

Acquire

Scheduling Scans: Set the Scan Frequency – Second Vessel

IC Add Vessel

Scan Schedule

Define the scan schedule for this vessel.

00 02 04 06 08 10 12 14 16 18 20 22 00

Add Scans to Schedule

☒ Create new schedule with scans at intervals of 2 Hours

☐ Add to existing schedule Next scan at 15:05, scanning every 2 hours

☐ Create new schedule with advanced scheduling options

☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

☒ Scan indefinitely

☐ Stop scanning 1 days, 0 hours after the first scan

☐ Stop scanning after 11/11/2020 14:00

For advanced scheduling options, double click the timeline on the main schedule page.

Back Next

Schedule Scans

- Create an independent second schedule if desired
- Add to an existing schedule to scan vessels sequentially at the same frequency
- Fine adjustments can be made by click and drag

Acquire

Scheduling Scans: Sampling Interval & Schedule


Sampling Interval : Do not over-sample (how fast is the process you want to quantify?)

Suggested Frequencies:

- 15-30 minutes interval for Phagocytosis or Immune Cell Killing
- 0.5-2 hours for Scratch Wound and Chemotaxis
- 2-4 hours for Proliferation and Cell Health
- 4-6 hours for Colony Formation and Spheroids

Acquire

Scheduling Scans: Advanced Scheduling



Advanced Scheduling

- Double clicking the timeline opens up advanced scheduling
- Allows manipulation of schedules


Acquire

Scheduling Scans: 2 Heat Threshold Warnings Maintain Cell Health with Optimal Temperature Conditions

IC Add Vessel

Scan Schedule

Define the scan schedule for this vessel.



Add Scans to Schedule

- ☐ Create new schedule with scans at intervals of 45 Minutes
- ☒ Add to existing schedule: Next scan at 19:45, scanning every 45 minu...
- ☐ Create new schedule with advanced scheduling options
- ☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

- ☒ Scan indefinitely
- ☐ Stop scanning: 1 days, 0 hours after the first scan
- ☐ Stop scanning after: 12/06/2020 20:00

Warning - For optimal system performance, there should be no more than 12 hours of scanning in a 24 hour period.

Reason - You are currently exceeding this recommendation.

How to Fix

- Reduce the scanning frequency for one or more scheduled vessels, or
- Reduce the number of vessels, or
- Reduce the number of locations imaged in a given vessel.

Back Next

- For every 1 minute of scanning, there should be at least 1 minute of non-scanning
- 12 hours of scanning per day
- Reduce scan frequency if scanning in excess of guidelines

Recommendations in Orange : Do NOT Ignore

Acquire

Scheduling Scans: 2 Heat Threshold Warnings

Maintain Cell Health with Optimal Temperature Conditions

Scan Schedule

Define the scan schedule for this vessel.

00 02 04 06 08 10 12 14 16 18 20 22 00

Add Scans to Schedule

- ☐ Create new schedule with scans at intervals of 3 Hours
- ☒ Add to existing schedule: Next scan at 02:00, scanning every 3 hours
- ☐ Create new schedule with advanced scheduling options
- ☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

- ☒ Scan indefinitely
- ☐ Stop scanning: 1 days, 0 hours after the first scan
- ☐ Stop scanning after: 03/07/2020 23:00

Warning - For optimal system performance, the total duration of any one scan group should not exceed 45 minutes.

Reason - You are currently exceeding this recommendation.

How to Fix - Split these scan groups into multiple scan groups.

Back Next

- Scan groups should not exceed 45 minutes
- Split long scan groups into two groups

Recommendations in Orange : Do NOT Ignore

Acquire


Scheduling Scans: 2 Scan Overlap/Duration Warnings

Avoid Conflicts in Scheduling

IC Add Vessel

Scan Schedule

Define the scan schedule for this vessel.



Add Scans to Schedule

☒ Create new schedule with scans at intervals of

☐ Add to existing schedule

☐ Create new schedule with advanced scheduling options

☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

☒ Scan indefinitely

☐ Stop scanning after the first scan

☐ Stop scanning after

Warning - The 12:55 AM scan may not run if the prior scan runs longer than shown above.

Reason - Scan durations are estimates. If the 12:26 AM scan overruns, the 12:55 AM scan may not run.

How to Fix - To reduce the chance of this occurring, increase the time between the 12:26 AM and 12:55 AM scans.

[Back](#) [Next](#)

- Scan durations are estimates, if a vessel overruns, the subsequent group may be skipped
- Placing scan groups too close together will lead to the Incucyte showing the scans as orange bars
- Drag orange bars to add more time in between scans

Recommendations in Orange : Do NOT Ignore

Acquire


Scheduling Scans: 2 Scan Overlap/Duration Warnings

Avoid Conflicts in Scheduling

IC Add Vessel

Scan Schedule

Define the scan schedule for this vessel.



Add Scans to Schedule

☒ Create new schedule with scans at intervals of

☐ Add to existing schedule

☐ Create new schedule with advanced scheduling options

☐ Reserve tray location and defer scheduling until later

Total Duration of Experiment

☒ Scan indefinitely

☐ Stop scanning after the first scan

☐ Stop scanning after


You have scan times that overlap. You must change the times of the scans so that no overlaps (shown in red in the timeline) occur.

[Back](#) [Next](#)

- Scans cannot be scheduled to occur at the same time
- The Incucyte will show overlapping scans as red bars
- Drag red bars to add more time in between scans
- Will not let you ignore

Acquire

Scheduling Scans: Drawer Layout Summary



- Vessel selected in 'blue'
- View Scan Properties/Notebook/Analysis
- Schedule is white for this vessel, other vessels are greyed out
- Add, Remove and Edit vessels using the icons above the tray layout

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion



Visualize


View Images

4 Steps to Image Analysis


Acquire Images

Scan Properties	
Vessel Type	24-well Corning Falcon
Scan Type	Standard
Image Channels	<input checked="" type="checkbox"/> Phase <input checked="" type="checkbox"/> Green Acquisition Time (ms) 300
Objective	20x
Scan Duration	2 min (estimated)
Number of Daily Scans	24


View Images



Process Images




Export Data




Incucyte® : A Guided Interface


Schedule Scans




Vessel View



Analysis Definition



Graphing and Exporting



Visualize

Search for Experiments by Label/user/etc.

Incucyte Toxicity Panel Assay (Green Reagents) - IncuCyte

Schedule View Manage Device Archive Status

Next Scan: 4 min

Enter text to search...

Analyses	Vessel Name	Owner	Last Scan	Scan Type	Vessel ID
▼	Toxicity Panel Assay (Green Reagents)	MeaganR	1/2/2017 1:13 PM	Standard	1013

Analysis Definition Name	Analysis Type	Creator	Date Completed	Analysis ID	Analysis Notes
CytoTox LO	Basic Analyzer	LindyO	1/3/2017 4:58 PM	4	
Cytotox Green	Basic Analyzer	MeaganR	1/3/2017 1:52 PM	3	
Caspase 3/7	Basic Analyzer	MeaganR	1/3/2017 1:41 PM	2	
Annexin Green	Basic Analyzer	MeaganR	1/3/2017 1:31 PM	1	

Plate Map: 12x8 grid (A-H, 1-12) showing well status.

Vessel Name: Toxicity Panel Assay (Green Reagents)

Scan Type: Standard

Cell Type: HT-1080 Lenti NLR

Image Channels: ☒ Phase ☒ Green ☒ Red

Magnification: 10x

Plate Map: [Yes](#)

Vessel Type: 96-well Corning


Vessel ID: 1013

Vessel Notes:

Visualize


Vessel (Image) Display

Packed View



Quickly assess trends in data in default 'Packed View'


Unpacked View



Visualize location of images in the well

Visualize

View and Interact with Images




Calibrated Images

- Compare images over time and between instruments
- Images made up of pixels with values used to define analysis
- Fluorescent pixel intensities are reported in calibrated units, GCU (Green), OCU (Orange) etc. Hover over any pixel to find these values.

Visualize

View and Interact with Images



Visualize images

- Adjust Brightness/Contrast for phase-contrast
- Adjust min/max intensity settings for fluorescence
- Adjusting these values will not change the pixel values of the image

Visualize

Autoscale


What is it?

- Sets a range of min/max intensities to account for changes in biological samples across different samples and time
- A 'halo' or 'cross' artifact might appear when imaging samples with dim signal

Minimum Green
intensity in the
Well/Vessel
= Black




Maximum Green
intensity in the
Well/Vessel



Visualize


Autoscale



Visualize

Vessel – Autoscale

Applies an autoscale across entire vessel based on min/max brightness across entire plate




Ideal for experiments with the same reagent.

In this experiment, 3 sections of the plate were treated with different reagents of different fluorescence intensity; therefore, vessel autoscale is not recommended

Visualize

Well – Autoscale

Applies an autoscale for each individual well based on min/max brightness in each individual well



Ideal for experiments with different reagents and to determine the baseline intensities of individual wells.


In this example, well autoscale is recommended

Visualize


Turn off Autoscale (recommended for movie making)

Use the default autoscale settings to find a base line by turning off autoscale. Adjust min and max intensity values for optimal visualization (does not alter pixel intensity)

Auto-scale applied



Manually adjusted min/max settings




Example of same image (1 with auto-scale, 1 with manual adjustment). Pixel intensity remains the same despite changing the visualization range.

Analogy: Think of playing the same song at different volumes.


Why is the Near-IR Channel Image Blue?

- The NIR emission wavelength for the Incucyte® SX5 is a deep red which does not provide visual contrast with the emission wavelength of the orange channel
- These are spectrally separated optically but the perceived colors are very close visually
- We use a blue pseudocolor to represent NIR to enable easy differentiation between channels, particularly when displayed as a multichannel overlay image



Spectral Unmixing

- Most fluorophores have broad and uneven excitation and emission spectra
- Bleedthrough is the consequence of the excitation of one fluorophore by more than one excitation channel, leading to an emission readout in multiple channels
- We use Spectral Unmixing to separate and remove the contribution of undesirable emission signal between fluorophores, providing an accurate representation of the contribution of each fluorophore, imaged in the same sample




<https://www.microscopyu.com/techniques/confocal/spectral-imaging-and-linear-unmixing>

Spectral Unmixing

- Fluorescent Dye Optimization Experiment
 - To test the extent of spectral overlap of a fluorophore
 - Prepare wells with individual reagents only (Red or Green)
 - Image all wells with both Red and Green channels
 - Check for signal in both channels, if present in the opposite channel, Spectral Unmixing is required
 - For example, for a Red fluorophore with signal appearing in the Green channel
 - Increase the %R contributes to G until no signal is seen in the Green channel


Incucyte® Spectral Unmixing Tool




Visualize

Spectral Unmixing Example


1. Red only fluorophore – imaged in both red and green channel




3. Increase the % R contributes to G until there is little/no green signal left



2. Turn off the Red channel and determine if there is any green signal




4. Adjust the min/max intensities to remove remaining green background



Visualize


Nuclight Orange Spectral Unmixing with the Incucyte® SX5

No Spectral Unmixing



- Orange and Green overlap makes the color appear yellow
- There is no green fluorescent protein expressed in this cell
- Green signal is coming from the Orange bleeding into Green


Spectral Unmixing applied



- After Spectral Unmixing is applied, the bleed through Green signal is removed and the overlapping image appears as its true Orange color

Spectral Unmixing Using the Green/Orange/NIR Optical Module

- Using 3 fluorescent channels gives 6 potential Spectral Unmixing values to enter - the software lists combinations from most likely to least likely to occur
 - Bleed through from Orange reagents into the Green channel and NIR reagents into the Orange channel are most common
 - Green bleeding into Orange, or Orange bleeding into NIR is not as common but does happen
 - Example: Caspase 3/7 Green slightly bleeds into the Orange channel and requires a value of 1.5% G contributes to O as shown
 - Due to increased spectral separation, bleed through between NIR and Green channels is unusual, but still possible



The screenshot shows the 'Spectral Unmixing' interface with six rows of input fields. Each row represents a contribution from one channel to another. The first two rows show 0% contribution, while the third row shows 1.5% contribution from Green to Orange.

Channel	Contributes to	Value (%)
% O	contributes to G	0
% N	contributes to O	0
% G	contributes to O	1.5
% O	contributes to N	0
% G	contributes to N	0
% N	contributes to G	0

Spectral Unmixing – Recommended Ranges for Incucyte® Reagents

Channel	Reagent	% O to G	% N to O	% G to O	% O to N	% G to N	% N to G
Green	Angiogenesis Prime Kit			-		-	
	Annexin V Green			-		-	
	Caspase 3/7 Green			1-2%		-	
	CytoLight Lenti			-		-	
	CytoLight Rapid Green			-		-	
	Cytotox Green			-		-	
	FabFluor IgG1, IgG2a, IgG2b			-		-	
	NucLight Green (BacMam?/lenti)			-		-	
	pHrodo e.coli, sa, zym Green			NR*		-	
Orange	Annexin V Orange	5-7%			-		
	NeuroBurst	-			-		
	NeuroLight	10-12%			-		
	NucLight Orange	10-12%			-		
	pHrodo cell labeling Orange	NR*			-		
	pHrodo e.coli, sa, zym Orange	NR*			-		
	MMP reagent	5-7%			-		
NIR	Annexin V NIR		1-2%				-
	NucLight NIR		1-2%				-
	NucLight Rapid NIR		1-2%				-
Green and Orange	FUCCI G/O	12-14%		-			

*Not Recommended for multiplexing

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion




Create an Analysis Definition

4 Steps to Image Analysis


Acquire Images

Scan Properties	
Vessel Type	24-well Corning Falcon
Scan Type	Standard
Image Channels	<input checked="" type="checkbox"/> Phase <input checked="" type="checkbox"/> Green Acquisition Time (ms) 300
Objective	20x
Scan Duration	2 min (estimated)
Number of Daily Scans	24


View Images



Process Images




Export Data




Incucyte® : A Guided Interface


Schedule Scans




Vessel View



Analysis Definition



Graphing and Exporting




Defining Analysis Masks

What is a mask?

- An area of interest (cells, nuclei, etc.) within an image
- Define by eye which pixels in an image correspond to the desired object by adjusting analysis parameters
- The same analysis is applied to all selected images

When do I make a new mask?


- When changing:
 - Cell type
 - Fluorophore
 - Magnification



The pseudo yellow color in this image indicates which phase areas of the image are of interest (cells) and should be quantified in the data analysis


Image Processing: 3 Steps to Image Analysis on the Incucyte® SX5

1. Select Representative Images




Include images from beginning/end, untreated/treated, dim/bright

2. Preview Analysis Definition




Adjust parameters to create mask

3. Select Time Points/Wells



Run analysis definition on all images/time points of choice

Image Processing: New or Existing Analysis



Analysis

- For a new experiment choose a New Analysis Definition
- If a repeat of a previous experiment use or copy an Existing Definition

Image Processing: Analysis Type




Choose the Appropriate Analysis Type

- The Basic Analyzer can be used to analyze assays such as proliferation, cell health and immune cell killing
- Cell-by-Cell can also be used for these assays yielding per cell information
- Certain imaging modes will predetermine the analysis type, e.g. Spheroid

Process

Image Processing: Select Image Channels




The screenshot shows a software window titled "Launch Analysis" with a sub-header "Image Channels". Below the header, a message states: "The Basic Analyzer requires at least one channel to be selected for analysis." Under the "Image Channel" section, four channels are listed with checkboxes: "Phase" (checked), "Green" (checked), "Orange" (checked), and "NIR" (checked). Under the "Overlap" section, there is a toggle switch labeled "No" and "Yes", with the "Yes" option selected. At the bottom of the window, there are "Back" and "Next" buttons, and a series of seven dots indicating the current step in the process.

Choose Appropriate Channels

- You do not have to analyze all channels, e.g. phase for ICK
- Overlap option quantifies fluorescent coincidence of objects positive for two or more fluorescent channels

Process

Image Processing: Select Channels for Coincident Object Analysis




The screenshot shows a software window titled "Launch Analysis" with a sub-header "Coincident Objects". Below the sub-header is the instruction: "Quantify the extent of overlap between fluorescent objects within each image." Under the "Overlap Mask" section, there are four checkboxes: "Green + Orange" (checked), "Orange + NIR" (unchecked), "Green + NIR" (checked), and "Green + Orange + NIR" (unchecked). At the bottom of the window, there is a "Back" button on the left, a series of seven dots in the center (the fourth dot is blue, indicating the current step), and a "Next" button on the right.

Quantify Fluorescent Overlap

- If desired select relevant channels for coincident object analysis
- Coincident object analysis is dependent on the initial analysis of each independent channel

Process

Image Processing: Choose an Image Set




Select Appropriate Images

- Select a range of images
- Should include the extremes of your experiment and in between

Image Processing: Selection of Representative Images

- Images representing the variability inside an experiment (cell morphology, density, fluorescence intensity, background)
- HD-phase
 - Different confluency level (10-80%range)
 - Different cell morphology (differentiation, spreading, death)
 - Images with debris/scratches if representative
- Brightfield (Spheroid and Organoid scan types)
 - Different size and density
 - Different background
- Fluorescence
 - Different intensity after correct setting of Min/Max
 - Different background

Image Processing: Preview the Current Image



Always Preview First

- The software uses the images to train a mask
- The user can then make manual adjustments

Image Processing: Adjust Masking

Launch Analysis

Analysis Definition

Refine the Analysis Definition settings using Preview to analyze the Image Set.

Metrics for Image D6-1 @ 2d 2h 0m <3 Color_Cell Cycle_Annexin NIR_Demo data set>

Image Channel	Confluence (%)	Count (Per Image)	Avg Area (µm²)	Avg Eccentricity	Avg Mean Inter...	Avg Integrated Int...	Total Integrated Int...
Confluence	57.968	367	3579.3	0.6956	N/A	N/A	N/A
Green Object	2.4482	425	130.54	0.6836	3.6564	584.29	2.4832E+05
Orange Object	0.6764	169	90.691	0.6872	4.2056	488.67	8.2585E+04
NIR Object	2.3118	881	59.463	0.6175	7.1580	789.02	6.9513E+05
Green + Orange	0.0797	32	56.436	N/A	N/A	N/A	N/A

Back Next

Manual Changes

- Go from top to bottom, changing parameters sequentially and checking
- Eccentricity is most often used when filtering clusters, not monolayers, due to variance in cell shape

Image Processing: Fluorescence Analysis

Launch Analysis

Analysis Definition

Refine the Analysis Definition settings using Preview to analyze the Image Set.

Analysis Definition

Phase

Green

Object Name

Segmentation

Top-Hat

Radius (µm) 20.000

Threshold (GCU) 2.0000

Edge Split On

Edge Sensitivity

Cleanup

Hole Fill (µm²) 0.0000

Orange

NIR

Green + Orange

Green + NIR

Whole Well

Layers

Image Channels

☒ Phase

☒ Green

☒ Orange

☒ NIR

Analysis Masks

Mode BLEND OVERLAY

Channel Masks

Metrics for Image D6-1 @ 2d 2h 0m <3 Color_Cell Cycle_Annexin NIR_Demo data set>


Image Channel	Confluence (%)	Count (Per Image)	Avg Area (µm²)	Avg Eccentricity	Avg Mean Intensity...	Avg Integrated Int...	Total Integrated Int...
Confluence	57.968	367	3579.3	0.6956	N/A	N/A	N/A
Green Object	1.9755	366	122.31	0.6777	3.8588	579.01	2.1192E+05
Orange Object	0.6589	168	88.872	0.6853	4.2156	482.34	8.1032E+04
NIR Object	2.2036	849	58.817	0.6133	7.3146	806.00	6.8429E+05
Green + Orange	0.0712	28	57.590	N/A	N/A	N/A	N/A

Back Next

Increased Parameters

- Always start with Top-Hat, selecting an appropriate radius
- Select a threshold e.g., GCU or OCU based on an unlabeled control

Image Processing: Top-Hat Background Subtraction for Fluorescence



Top-Hat:

- Radius larger than that of the largest object of interest
- Separation into 2 components: background and real signal
- Threshold applied to include only real signal in the mask

How does it work:






Image Processing: Fluorescence Analysis

Launch Analysis

Analysis Definition

Refine the Analysis Definition settings using Preview to analyze the Image Set.

Analysis Definition

Phase

Green

Object Name

Segmentation

Top-Hat

Radius (µm) 20.000

Threshold (GCU) 2.0000

Edge Split On

Edge Sensitivity

Cleanup

Hole Fill (µm²) 0.0000

Orange

NIR

Green + Orange

Green + NIR

Whole Well

Layers

Image Channels

☒ Phase

☒ Green

☒ Orange

☒ NIR

Analysis Masks

Mode BLEND OVERLAY

Channel Masks

Metrics for Image D6-1 @ 2d 2h 0m <3 Color_Cell Cycle_Annexin NIR_Demo data set>

Image Channel	Confluence (%)	Count (Per Image)	Avg Area (µm²)	Avg Eccentricity	Avg Mean Intensity...	Avg Integrated Int...	Total Integrated Int...
Confluence	57.968	367	3579.3	0.6956	N/A	N/A	N/A
Green Object	1.9755	366	122.31	0.6777	3.8588	579.01	2.1192E+05
Orange Object	0.6589	168	88.872	0.6853	4.2156	482.34	8.1032E+04
NIR Object	2.2036	849	58.817	0.6133	7.3146	806.00	6.8429E+05
Green + Orange	0.0712	28	57.590	N/A	N/A	N/A	N/A

Back Next

Analyze each channel

- For nuclear labels keep Edge Split on
- For membrane and cytoplasmic labels turn Edge Split off
- For nuclear counts a size filter is essential to exclude vesicles and debris

Process

Image Processing: Select Wells and Times to Analyse


The screenshot shows a software window titled 'Launch Analysis' with a sub-header 'Scan Times and Wells'. Below the header is a instruction: 'Select the scan times and images to analyze along with the option to analyze future scans.' The interface is divided into two main sections: 'Select Scans' on the left and 'Select Wells' on the right. In the 'Select Scans' section, there is a list of time points from '0d 0h 0m' to '2d 2h 0m' in 2-hour increments. The first six time points (0d 0h 0m to 0d 22h 0m) are highlighted in blue, and there is an unchecked checkbox for 'Analyze Future Scans'. In the 'Select Wells' section, there is a grid of 96 wells (8 rows labeled A-H and 12 columns labeled 1-12). Each well is represented by a circle with a small square in the center. Wells A1 through A6, B1 through B6, C1 through C6, D1 through D6, E1 through E6, F1 through F6, G1 through G6, and H1 through H6 are all selected, indicated by blue circles. Wells A7 through A12, B7 through B12, C7 through C12, D7 through D12, E7 through E12, F7 through F12, G7 through G12, and H7 through H12 are not selected, indicated by white circles. At the bottom of the window, there are 'Back' and 'Next' buttons, and a progress indicator showing 8 steps with the 4th step (corresponding to 'Select Wells') being the current active step.

Choose and Include Appropriate Wells

- Analyse wells of interest
- Exclude individual timepoints if necessary

Process

Image Processing: Name the Analysis



The screenshot shows a software window titled "Launch Analysis". Inside, the main heading is "Save and Apply Analysis Definition". Below this, a instruction reads: "Assign the Analysis Definition a name and, if desired, add notes." There are two input fields: "Definition Name" and "Analysis Notes". The "Definition Name" field contains the text "HT-1080 Cell Cycle Green Orange Annexin V NIR". The "Analysis Notes" field is empty. At the bottom of the window, there is a "Back" button on the left, a "Next" button on the right, and a series of seven dots in the center, with the fourth dot from the left being highlighted in blue, indicating the current step in the process.

Descriptive Name

- Try to be clear with naming and contain details about the analysis
- Possible to share analyses between users

Process

Image Processing: Summary Screen

Launch Analysis

Summary

Verify that the information below is correct before launching the Analysis.

Analysis Type	
Type	Basic Analyzer
Image Channels	
Channels	Phase, Green, Orange, NIR
Coincident Objects	
Overlap Masks	Green + Orange, Green + NIR
Image Set Selection	
Images Previewed	2
Select Times and Wells	
Analyze Future Scans	No
Wells to Analyze	48
Scan Times to Analyze	14
Analysis	
Definiton Name	HT-1080 Cell Cycle Green Orange Annexin V NIR
Notes	
Vessel Information	
ID	481
Name	3 Color_Cell Cycle_Annexin NIR_Demo data set
Vessel Type	96-well TPP
Scan Type	Adherent Cell-by-Cell

[Back](#) [Finish](#)

Check the Analysis Looks Correct

- If changes need to be made, click the tab and it will link back
- If all is correct click Finish and the analysis will start

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion




Export Data and Images

4 Steps to Image Analysis


Acquire Images

Scan Properties	
Vessel Type	24-well Corning Falcon
Scan Type	Standard
Image Channels	<input checked="" type="checkbox"/> Phase <input checked="" type="checkbox"/> Green
	Acquisition Time (ms) 300
Objective	20x
Scan Duration	2 min (estimated)
Number of Daily Scans	24


View Images



Process Images




Export Data




Incucyte® : A Guided Interface


Schedule Scans




Vessel View



Analysis Definition



Graphing and Exporting



Export: Search Analyses

- Analyses associated with a vessel are indicated by a blue arrow
- Analyses saved by Analysis Definition Name
- Right click on Analysis Definition to view details
- Directly go to graphing window
- Analysis Notes specify details of analysis definition
- Details in the right pane show wells/image channels analyzed

The screenshot displays the Incucyte software interface. At the top, there is a navigation bar with icons for Schedule, View, Manage, Device, Archive, and Status. A search bar is located below the navigation bar. The main area is divided into two panes. The left pane shows a table of analyses, and the right pane shows details for the selected analysis.

Analyses Table:

Analyses	Vessel Name	Owner	Last Scan	Scan Type	Vessel ID
	Toxicity Panel Assay (Green Reagents)	MeaganR	1/2/2017 1:13 PM	Standard	1013

Analysis Definition Table:

Analysis Definition Name	Analysis Type	Creator	Date Completed	Analysis ID	Analysis Notes
CytoTox LO	Basic Analyzer	LindyO	1/3/2017 4:58 PM	4	
Cytotox Green	Basic Analyzer	MeaganR	1/3/2017 1:52 PM		
Caspase 3/7	Basic Analyzer	MeaganR	1/3/2017 1:41 PM		
Annexin Green	Basic Analyzer	MeaganR	1/3/2017 1:31 PM		

Right Pane Details:




- Analysis Type:** Basic Analyzer
- Analysis Definition:** Caspase 3/7
- Analysis ID:** 2
- Vessel Name:** Toxicity Panel Assay (Green Reagents)
- Scan Type:** Standard
- Cell Type:** HT-1080 Lenti NLR
- Image Channels:** ☒ Phase, ☒ Green, ☒ Red
- Magnification:** 10x
- Plate Map:** [Yes](#)
- Vessel Type:** 96-well Corning
- Vessel ID:** 1013
- Analysis Notes:** (highlighted)

Well Plate Map:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Data Export



Export: Graph – View/Export



- Select time points, wells, metric and grouping of interest
- Possible to build custom metrics, including data normalisation

Data Export


Export: Graph – View/Export



- Customize graphs to adjust colors, legends, axes etc.
- Combine graphs and compare datasets using Drag and Drop

Data Export

Export: Graph – View/Export



The Graphing Export dialog box provides options for exporting the data. The 'Layout' section has three radio buttons: 'Show each scan as a single row in one large table.', 'Column by column: A1, B1, C1, ... A2, B2, C2, ...', and 'Row by row: A1, A2, A3, ... B1, B2, B3, ...'. The 'Destination' section has three radio buttons: 'Clipboard', 'All scans in one file' (selected), and 'Each scan in a separate file'. The 'Other Options' section includes checkboxes for 'Include experiment details in header', 'Include', and 'Break data down into individual images'. There are also options for 'Fill holes in the data with the following characters' and 'Calculate Error Per Image' or 'Calculate Error Per Well'. A 'Grouping' dropdown menu is set to 'None'. A yellow box highlights the 'Grouping' dropdown menu, and a yellow line points to it from the 'None' option in the 'Microplate Graph' window.



- Export Microplate Graphs to assess data trends and outliers
- Export data for 3rd party statistical analysis, select grouping as 'None'

Export: Images and Movies

- Export “As Displayed”
 - Presentations/Publications
 - Adjust Brightness/Contrast for phase-contrast
 - Adjust min/max intensity settings for fluorescence (remove autoscale)
- Export “As Stored”
 - Analyze images for 3rd Party analysis packages
 - Phase: 8-bit raw
 - Fluorescence: Uncalibrated raw 16-bit or calibrated 32-bit floating point



Data Export

Export: Images and Movies



- Enter Guided Interface and select Export Type
- Select images of interest, channels and masks and adjust fluorescence scaling

Export: Images and Movies



- Select Single Movie or Series of Images and select timepoints
- Adjust scaling, frame rate, timestamp and legend

Agenda

Introduction and Applications

Incucyte® Hardware and Best Practices

Software Overview

Acquire

Visualize

Process

Data Export

Conclusion



Conclusion

Resources – Publications Library



Over 3900 peer reviewed publications and growing!

Explore the Incucyte® publications library - an up-to-date and fully searchable resource!

[Incucyte® Publication Library](#)

Conclusion

Incucyte® Support



Mathias Lucas

FAS, France

mathias.lucas@sartorius.com

+33 6 42 02 62 28

Additional Resources:

- Quotes, Orders:

Olivier Bruno

olivier.bruno@sartorius.com

+33 33 6 60 27 46 53

- Assays, Reagents, Protocols, Software:

www.essenbioscience.com

askascientist@sartorius.com

- Technical Support and Hardware:

EUincucyte.support@sartorius.com

Thank You

Mathias Lucas
Field Application Scientist
Mobile +33 6 42 02 62 28
mathias.lucas@sartorius.com

The Sartorius logo is displayed in a bold, black, sans-serif font. The letters are closely spaced, and the 'S' and 'A' are particularly prominent. The logo is centered within a bright yellow rectangular area that occupies the bottom right portion of the slide.



Simplifying Progress

Visit
sartorius.com/incucyte

SARTORIUS

Main sources of pipetting errors according to ISO 8655-2

Greatest impact on the pipetting results - up to 50%:

- Leaky piston and cylinder system of the pipette
- Poor tip fitting (tip is leaking)

Medium impact to the pipetting results - up to 4%:

- Re-use of the tip
- Failure to wipe the tip against the vessel wall (10-15mm upwards)
- Unstable humidity of the pipetting environment
- No pre-rinsing of the tip
- Uneven rhythm and timing of pipetting
- Immersion depth of the tip (2-3mm) and angle of the pipette during pipetting

Minor impact on the pipetting results, up to 0,5%:

- Inconsistent piston movement
- Difference in temperature between the pipette, tip, liquid and room temperature (up to 0,3%/C°)